SENO MARK ULTRACOR™ MRI

SENO MARK ULTRACOR™ MRI 14 Gauge Tissue Marker And Introducer Kit

ENGLISH

INSTRUCTIONS FOR USE

REF	SMUCMRI14GT (Titanium Marker)
	SMUCMRI14GSS (Stainless Steel Marker)
	UCMRINTLOC (Introducer Kit – packaged separately)
	14 Gauge Tissue Marker and Introducer Kit

Caution: Federal (U.S.A.) law restricts these devices to sale by or on the order of a physician.

DESCRIPTION

The SenoRx Inc. SENOMARK ULTRACOR™ MRI Tissue Marker consists of a disposable 14 gauge beveled needle Applicator containing:

- 3 resorbable polyglycolic acid (PGA) pads, with a wireform embedded in the middle pad;
- polyethylene glycol (PEG) plug at the distal end.

The wireform is intended for long-term radiographic marking of the tissue site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The Applicator has a beveled 13 cm/14G needle with 0.5 cm depth marks and a locking plunger. The pads are deployed from the beveled needle tip into the tissue site.

The ULTRACOR® MRI Introducer Kit consists of:

- a titanium Trocar with a plastic Sheath;
- a plastic Obturator;
- a Needle Guide Block.

The titanium Trocar and plastic Sheath are used to access the tissue site. The Needle Guide Block is compatible with Invivo¹ Grid Localization System.

INTENDED USE

The SENOMARK ULTRACOR™ MRI is intended to radiographically and sonographically mark breast tissue prior to or during a breast biopsy procedure.

CONTRAINDICATIONS

This device is not intended for use except as indicated.

WARNINGS

- **Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.**
- **The SENOMARK ULTRACOR™ MRI is not recommended for use in patients with breast implants.**
- **Do not use in the presence of infection.**

PRECAUTIONS

- Store at temperatures below 25°C (77°F).
- These devices should only be used by physicians trained in percutaneous breast biopsy procedures.
- Do not use these products if the sterile barrier has been previously opened or if the package is damaged.
- These products are provided sterile, and are intended for single use only. **DO NOT RESTERILIZE.**
- Do not use if needle is bent and/or tip is damaged.
- Verify Sheath does not move when removing Trocar/Obturator.
- Verify Sheath does not move when placing Applicator within the Sheath.
- Ensure that all pads are dispensed.
- Published studies with comparably sized 316L stainless steel biopsy markers have shown no marker movement and

insignificant marker heating when tested in an MRI System with a 1.5T field strength.

COMPLICATIONS

Potential complications (e.g. infection) that may be associated with the use of the SENOMARK ULTRACOR™ MRI and ULTRACOR® MRI Introducer Kit are similar to those associated with the use of other tissue marking devices.

HOW SUPPLIED

The SENOMARK ULTRACOR™ MRI and ULTRACOR® MRI Introducer Kit are provided sterile, and are intended for single use only.

DIRECTIONS FOR USE (Refer to Figure 1)

ULTRACOR® MRI Introducer Kit

1. Inspect the packaging to insure that package integrity has not been compromised. The product is sterile unless the seal is broken.
2. Using standard aseptic technique, open the package and remove the Titanium Trocar and Sheath (N), Needle Guide Block (J) and Obturator (M).
3. Identify the target lesion or site in the breast.
4. Use the Needle Guide Block with the MRI Devices Grid. Make a skin nick.
5. Remove the Tip Protector (K) and remove adapter (I) from the Trocar and Sheath assembly and position the stop (L) on the Sheath. Using the 0.5 cm depth markings, insert the Trocar and Sheath assembly into the breast, directing the Trocar tip to the target site. Trocar may be replaced with Obturator for imaging.
6. Re-image the breast to verify placement of the Trocar tip at the target site. Modify position, if needed.


SENO MARK ULTRACOR™ MRI Applicator placement using the ULTRACOR® MRI Introducer Kit

1. Stabilize the Sheath while removing the Trocar/Obturator and inserting the Applicator (B) into the Sheath. Utilize the depth markers and ULTRACOR® MRI logo (F) to confirm depth and orient needle bevel.


SENO MARK ULTRACOR™ MRI Applicator placement using the ENCOR® MRI Introducer Kit

Note: refer to ENCOR® MRI Introducer Kit ECMRINTLOC instructions for use and description of components.

1. Remove the Obturator/Probe from the cannula.
2. Place the Adapter onto the needle of the applicator as shown in the figure below. Placing it firmly against the handle of the applicator.

3. Insert the Applicator into the cannula until the front portion enters the seal of the cannula and stops at the first edge.

4. Unlock plunger by rotating clockwise 90°.
5. Deploy the marker by advancing plunger until it contacts barrel of applicator.
6. Remove Applicator and cannula simultaneously and place sterile gauze beneath the guide block to temporarily hold pressure during imaging.
7. Re-image the breast to verify placement of the wireform at the target site.
8. Hold compression as described by standard medical guidelines.

Pellet Deployment from SENOMARK ULTRACOR™ MRI Applicator



Note that the ULTRACOR® MRI logo (F) is aligned with the needle bevel.

1. Release the Plunger safety lock (G) by turning the Plunger (A) clockwise 90 degrees until the lock symbol is vertically aligned (H) with the logo.
2. Immediately deploy all of the contents of the Applicator by advancing the Plunger until it contacts the barrel of the Applicator, while maintaining the Applicator stable in position. Remove Applicator from Sheath and replace with Trocar or Obturator.
3. Re-image the breast to verify placement of the wireform at the target site.
4. Remove the Sheath and Trocar/Obturator, maintaining compression to the needle track to minimize bleeding. Hold compression until bleeding has ceased.
5. Dispose of the Applicator, Sheath, Obturator, Needle Guide Block and Trocar properly.


Figure 1.

SENOMARK ULTRACOR™ MRI Tissue Marker and Introducer Kit

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ENGLISH

This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

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© 2012 by C. R. Bard Inc. All rights reserved.
SenoRx, Inc.
Bard Limited